



Tummy tucking
is easy
 with

CRAYO

SPECIFICATION

Model	CoolLipo Super
The material of Crayo handle	Medical Silicone
Cryolipolysis temperature	5°C to 11°C adjustable
Power	900W
Vacuum	650mmhg
Heat temperature	37-45°C
Heat temperature time	3 minutes
Heat temp Advantages	Avoid the frostbit during Crayo treatment
Light & Wavelength	Red (630nm) Green (570nm) Blue (480nm)
Cavitation frequency	40Khz
RF Frequency	10Mhz
Lipo Laser Big Power	100mw/pcs
Lipo Laser quantity	42pcs
Lipo laser big power pads	Big pads*4 + small pads*2
Machine Display	8.4 inch touch screen
200 size Crayo handle	1pcs
150 size Crayo handle	1pcs
Cooling System	Semiconductor + water + air
Energy controlled	Yes
Power supply	110V-120V/220V-240V, 50/60Hz

